

IATROBEADS

(packing material for Liquid Chromatography)

“IATROBEADS” consists of porous fine silica gel uniformly in sphere and particle size to provide excellent separation performance.

Characteristics

	IATROBEADS 6RS-8005	IATROBEADS 6RS-8010	IATROBEADS 6RS-8060	IATROBEADS 6RS-80100
Apparent density (g/ml)	0.42	0.46	0.58	0.51
pH	7.1	7.1	7.1	7.1
Micropore volume (ml/g)	0.87	0.88	0.80	0.81
Average pore size (Å)	80	80	80	80
Average particle size (μ m)	5	10	60	100
Water content (Less then %)	0.15	0.3	1.5	1.5

I A T R O B E A D S
(Packing material for Liquid Chromatography)

" IATROBEADS " consists of porous fine silica gel uniformly in sphere and particale size to provide excellent separation performance.

I A T R O B E A D S 6 R S - 8 0 6 0



C H A R A C T E R I S T I C S

Apparent density	0.58 g/mL
pH	6.8
Water-soluble substances	0.15 %
Surface area	420 m ² /g
Micropore volume	0.85 mL/g
Mean diameter of micropores	80 A
Water content	Less then 2.0 %
Average particle size	60 μm

S e p a r a t i o n b y " I A T R O B E A D S "

Elution patterns of human erythrocyte glycolipids from IATROBEADS GRS-8060 Column.

Thin-layer chromatogram of glycolipids isolated by IATROBEADS GRS-8060 Column.



Solvent system:

CHCl₃:CH₃OH:H₂O

83:16:0.5 to

55:42:3 (linear gradient)

Peak 1. CDH- I 2. CDH- II 3. CTH- I
4. CTH- II 5. Glob- I

CHCl₃:CH₃OH:3.5M NH₄OH

65: 35 : 8


1. CDH, CTH Glob- I
2. CDH- I 3. CDH- II
4. CTH- I 5. CTH- II
6. Glob- I

Reference: Ando, S., Isobe, M. and Nagai, Y.
(1976) Biochim. Biophys. Acta 424, 98-105

Silica-gel (Iatrobeds) column chromatography

【 Experiment】 Silica-gel column chromatography

- (1) Suspend 80 gr. of Iatrobeds (6RS-8060) in C-M-W (80:15:0.5) , then agitate the suspension in ultrasonic water bath.
→ Air bubbles will continue being generated for long time. But, only one (1) to two (2) minute-agittation will enable to prevent from air bubbling in the column during experiment.
- (2) Pack the suspended Iatrobeds (solid volume : 160 ml) into a column (1.8 cm ϕ \times 60 cml), and hence rinse the Iatrobeds in it with 200 ml of the same solvent.
- (3) Dissolve 0.5 g of glycolipids mixture in 5 ml of same solvent and put it on the column.
- (4) Elute the glycolipids mixture by means of concentration gradient method using 500 ml of C-M-W (85:15:0.5) as the first solvent and 600 ml of C-M-W (20:80:5) as the second solvent.
- (5) Almost pure standard component from monohexylceramide to Forssmann glycolipid can be obtained component by component in onece column chromatography by analyzing each 5 ml of elute on thin layer chromatography one by one which brings the separation pattern shown below.




A : monohexyleceramide B : dihexyleceramide C : trihexyleceramide
D : globoside E : Forssmann glycolipid F : neutral lipids with more than 6 sugars

Fig. Separation of glycolipids by concentration gradient method

Reference: Dr. S. Ando – Tokyo Metropolitan Institute of Gerontology

Separation of Lipids

Column: Packed IATROBEADS 6RSP-8010 (4.6 × 250mm)



Sample: Lipids Mixture

Eluent: n-Hexane: Isopropyl Alcohol: 0.2M Acetic acid Buffer (pH 4.2)

8 : 8 : 1

Detector: 206nm

Flow rate: 2 ml/min

Equipment: Hitach model L-6000

Regeneration of Iatrobeads (6RS-8060)

The used Iatrobeads are regenerated by the following treatment procedure, and then can be used several times as almost new Iatrobeads.

- (1) Agitate the 0.1N-HCl suspension of Iatrobeads in the Erlenmeyer flask.
- (2) Filtrate the suspension on the funnel to remove the HCl solution.
- (3) Rinse the Iatrobeads on the funnel with distilled water.
- (4) Repeat the steps from (1) to (3) three times.
- (5) Then, add distilled water on the Iatrobeads in the funnel repeatedly until outflow water becomes neutral by checking with PH-test paper.
- (6) Activate the Iatrobeads treated for 30 minutes at 120°C in the drying oven.